CAT Critically Appraised Topic

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SEARCH TERMS

- 1) MeSH Database (PubMed): MeSH term: "surgical procedures, minimally invasive", "parathyroidectomy", "parathyroid hormone", "hyperparathyroidism, primary", "protein stability", "immunoassay", "point-of-care systems"
- 2) UpToDate Online version 22.2 (2014): "parathyroidectomy", "PTH assay", "primary hyperparathyroidism"

RELEVANT EVIDENCE/REFERENCES

Guidelines and Recommendations

NICHOLS JH, CHRISTENSON RH, CLARKE W, et al. Executive summary: The National Academy of Clinical Biochemistry Laboratory Medicine Practice Guideline: Evidence-based practice for point-of-care testing. Clinica Chimica Acta; 2007. 379(1-2): 14-28. **(8)**

NATIONAL KIDNEY FOUNDATION. K/DOQI clinical practice guidelines for bone metabolism and disease in chronic kidney disease. Am J Kidney Dis 2003; 42(4 Suppl 3): S1-201. **(12)**

Reviews

FRASER WD. Hyperparathyroidism. Lancet 2009; 374: 145-158. **(1)**

COLE DEC, WEBB S, CHAN P. Update on parathyroid hormone: new tests and new challenges for external quality assessment. Clinical Biochemistry 2007; 40: 585-590. **(4)**

SOUBERBIELLE JC, FRIEDLANDER G, CORMIER C. Practical considerations in PTH testing. Clinica Chimica Acta 2006; 366: 81-89. **(5)**

SOKOLL LJ, WIANS FH JR, REMALEY AT. Rapid intraoperative immunoassay of parathyroid hormone and other hormones: a new paradigm for point-of-care testing. Clinical Chemistry 2004; 50(7): 1126-1135. **(6)**

BARCZYNSKI M, KONTUREK A, HUBALEWSKA-DYDEJCZYK A, et al. Evaluation of Halle, Miami, Rome, and Vienna intraoperative iPTH assay criteria in guiding minimally invasive parathyroidectomy. *Langenbecks Arch Surg* 2009; 394: 843-849. **(7)**

CHIU B, STURGEON C, ANGELOS P. Which intraoperative parathyroid hormone assay criterion best predicts operative success? A study of 352 consecutive patients. *Arch Surg* 2006; 141: 483-488. **(10)**

LEIKER AJ, YEN TWF, EASTWOOD DC, et al. Factors that influence parathyroid hormone half-life: determining if new intraoperative criteria are needed. *JAMA Surg* 2013; 148(7): 602-606. **(13)**

HELLER KS, BLUMBERG SN. Relation of final intraoperative parathyroid hormone level and outcome following parathyroidectomy. *Arch Otolaryngol Head Neck Surg* 2009; 135(11): 1103-1107. **(14)**

SCHNEIDER DF, OJOMO KA, MAZEH H, et al. Significance of rebounding parathyroid hormone levels during parathyroidectomy. *JSR* 2013; 184: 265-268. **(15)**

REIHER AE, SCHAEFER S, CHEN H, et al. Does the final intraoperative PTH level really have to fall into the normal range to signify cure? *Ann Surg Oncol* 2012; 19: 1862-1867. **(16)**

Original Articles

GREENE AB, BUTLER RS, MCINTYRE S, et al. National trends in parathyroid surgery from 1998 to 2008: a decade of change. *J Am Coll Surg* 2009; 209(3): 332-343. **(3)**

IRVIN GL 3RD, DEMBROW VD, PRUDHOMME DL. Clinical usefulness of an intraoperative "quick parathyroid hormone" assay. *Surgery* 1993; 114(6): 1019-1022. **(9)**

IRVIN GL 3RD, SOLORZANO CC, CARNEIRO DM. Quick intraoperative parathyroid hormone assay: surgical adjunct to allow limited parathyroidectomy, improve success rate, and predict outcome. *World J Surg* 2004; 28(12): 1287-1292. **(10)**

RICHARDS ML, THOMPSON GB, FARLEY DR, et al. An optimal algorithm for intraoperative parathyroid hormone monitoring. *Arch Surg* 2011; 146(3): 280-285. **(11)**

WESTERDAHL J, BERGENFELZ A. Parathyroid surgical failures with sufficient decline of intraoperative parathyroid hormone levels: unobserved multiple endocrine neoplasia as an explanation. *Arch Surg* 2006; 141: 589-594. **(17)**

ENGLISH E, MCFARLANE I, TAYLOR KP, et al. The effect of potassium EDTA on the stability of parathyroid hormone in whole blood. *Ann Clin Biochem* 2007; 44: 297-299. **(18)**

TEAL TK, REED M, STEVENS PE, et al. Stability of parathyroid hormone *ex vivo* in haemodialysis patients. *Ann Clin Biochem* 2003; 40: 191-193. **(19)**

HOLMES DT, LEVIN A, FORER B, et al. Preanalytical influences on DPC IMMULITE 2000 Intact PTH assays of plasma and serum from dialysis patients. *Clinical Chemistry* 2005; 51(5): 915-917. **(20)**

ANDERSON NR, NICHOLAS J, HOLLAND MR, et al. Effect of a protease inhibitor on *in vitro* stability of intact parathyroid hormone. *Ann Clin Biochem* 2003; 40: 188-190. **(21)**

CAVALIER E, DELANAYE P, CARLISI A, et al. Stability of intact parathyroid hormone in samples from hemodialysis patients. *Kidney International* 2007; 72(3): 370-372. **(22)**

TERRIS DJ, WEINBERGER PM, FARRAG T, et al. Restoring point-of-care testing during parathyroidectomy with a newer parathyroid hormone assay. *Otolaryngology-Head and Neck Surgery* 2011; 145(4): 557-560. **(24)**

O'CONNELL DA, SEIKALY H, HARRIS JR. Central laboratory versus point of care testing in intraoperative monitoring of parathyroid hormone levels: cost comparison. *J Otolaryngol Head Neck Surg* 2008; 37(1): 91-97. **(25)**

CARNEIRO DM, IRVIN GL 3RD. New point-of-care intraoperative parathyroid hormone assay for intraoperative guidance in parathyroidectomy. *World J Surg* 2002; 26: 1074-1077. **(26)**

JARRIGE V, NIEUWENHUIS JH, VAN SON JPHF, et al. A fast intraoperative PTH point-of-care assay on the Philips handheld magnotech system. *Langenbecks Arch Surg* 2011; 396: 337-343. **(27)**

JACOBSON SR, VAN HEERDEN JA, FARLEY DR, et al. Focused cervical exploration for primary hyperparathyroidism without intraoperative parathyroid hormone monitoring or use of the gamma probe. *World J Surg* 2004; 28(11): 1127-1131. **(28)**

GAWANDE AA, MONCHIK JM, ABBRUZZESE TA, et al. Reassessment of parathyroid hormone monitoring during parathyroidectomy for primary hyperparathyroidism after 2 preoperative localization studies. *Arch Surg* 2006; 141(4): 381-384. **(29)**

STALBERG P, SIDHU S, SYWAK M, et al. Intraoperative parathyroid hormone measurement during minimally invasive parathyroidectomy: does it "value-add" to decision-making? *J Am Coll Surg* 2006; 203: 1-6. **(30)**

MERLINO JI, KO K, MINOTTI A, et al. The false negative technetium-99m-sestamibi scan in patients with primary hyperparathyroidism: correlation with clinical factors and operative findings. *Am Surg* 2003; 69(3): 225-229. **(31)**

GOLDSTEIN RE, CARTER WM, FLEMING M, et al. Unilateral cervical surgical exploration aided by intraoperative parathyroid hormone monitoring in patients with primary hyperparathyroidism and equivocal sestamibi scan results. *Arch Surg* 2006; 141: 552-559. **(32)**

BARCZYNSKI M, KONTUREK A, HUBALEWSKA-DYDEJCZYK A, et al. Utility of intraoperative bilateral internal jugular venous sampling with rapid parathyroid hormone testing in guiding patients with a negative sestamibi scan for minimally invasive parathyroidectomy-a randomized controlled trial. *Arch Surg* 2009; 394: 827-835. **(33)**

BARCZYNSKI M, KONTUREK A, CICHON S, et al. Intraoperative parathyroid hormone assay improves outcomes of minimal invasive parathyroidectomy mainly in patients with a presumed solitary parathyroid adenoma and missing concordance of preoperative imaging. *Clinical Endocrinology* 2007; 66: 878-885. **(34)**

SAKIMURA C, MINAMI S, HAYASHIDA N, et al. Can the use of intraoperative intact parathyroid hormone monitoring be abandoned in patients with hyperparathyroidism? *The American Journal of Surgery* 2013; 206(4): 574-577. **(35)**

TWIGT BA, VAN DALEN T, VOLLEBREGT AM, et al. The additional value of intraoperative parathyroid hormone assessment is marginal in patients with nonfamilial primary hyperparathyroidism: a prospective cohort study. *The American Journal of Surgery* 2012; 204: 1-6. **(36)**

Reference Works, Handbooks and Databases

BURTIS CA, ASHWOOD ER, BRUNS DE. *Tietz Fundamentals of Clinical Chemistry*, 6th edition. Saunders (USA), 2007: 721-723. **(2)**

Posters, "grey literature", presentations

SCHARNHORST V, VALKENBURG J, VOSTERS C, et al. Influence of preanalytical factors on the Immulite intact parathyroid hormone assay. *Clinical Chemistry* 2004; 50(5): 974-975. **(23)**

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	Assay			
	Quick-IntraOperative Intact PTH	Turbo Intact PTH	STAT-10-I-PTH	Quick-IntraOperative Bio-Intact PTH(1-84)
Manufacturer	Nichols Institute Diagnostics	Diagnostic Products Corporation	Future Diagnostics	Nichols Institute Diagnostics
Platform	Single-well luminometer	IMMULITE; IMMULITE 1000/2000	Strip Luminometer	Advantage
Automated system	No	Yes	No	Yes
Assay design				
Capture antibody	Goat polyclonal anti-PTH(39-84)	Goat polyclonal anti-PTH(44-84)	Goat polyclonal anti-PTH(39-84)	Goat polyclonal anti-PTH(39-84)
Detector antibody	Goat polyclonal anti-PTH(1-34)	Goat polyclonal anti-PTH(1-34)	Goat polyclonal anti-PTH(1-34)	Goat polyclonal anti-PTH(1-5)
Form of PTH measured	Intact	Intact	Intact	Bio-Intact
Label	Acridinium ester	Alkaline phosphatase	Isoluminol	Acridinium ester
Signal	Chemiluminescence	Chemiluminescence	Chemiluminescence	Chemiluminescence
Solid phase	Bead	Bead	Microtiter strips with Accusphere	Magnetic particles (streptavidin coated)
Single-patient packaging	Yes	No	No	No
Specimen volume, µL	200	100	100	250
Time to result (including specimen processing), min	~15	~15	~10	~12
Incubation conditions	Shake at 400 rpm at 45 °C	Intermittent agitation at 37 °C	Shake at 1400 rpm at room temperature	37 °C
Measuring range, ^a ng/L	6-1250	4-2500	6-2500	1.2-5000
Calibration	Each assay	Readjust every 2 weeks	Each assay	Readjust every 7 days
Intraassay CV (concentration) ^a	9.4% (29 ng/L)	6.7% (49 ng/L)	10.7% (41 ng/L)	Readjust after 7 days (same reagent set) or 1 month (same reagent lot)
Reference interval, ^a ng/L	8.0% (162 ng/L)	6.5% (247 ng/L)	4.9% (255 ng/L)	5.4% (30 ng/L)
	10-65	12-72	10-65	4.0% (271 ng/L)
				15-65

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
English

Intended use

Immunoassay for the in vitro quantitative determination of intact parathyroid hormone in human serum and plasma for the differential diagnosis of hypercalcemia and hypocalcemia. This assay can be used intraoperatively. The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and cobas e immunoassay analyzers.

Summary

Parathyroid hormone (PTH) is formed in the parathyroid glands and secreted into the blood stream. Intact PTH consists of a single polypeptide chain containing 84 amino acids and has a molecular weight of approx. 9500 daltons.

The biologically active N-terminal fragment has a half-life of only a few minutes. Selective measurement of the (mainly) intact parathyroid hormone permits direct ascertainment of the secretory activity of the parathyroid glands.^{1,2}

PTH, together with vitamin D and calcitonin, brings about mobilization of calcium and phosphate from the skeletal system and increases the uptake of calcium in the intestine and the excretion of phosphate via the kidneys. The constancy of the blood calcium level is ensured by the interaction of PTH and calcitonin. The secretion of PTH is inhibited by high calcium concentrations and promoted by low calcium concentrations.

Parathyroid gland disorders lead to elevated or depressed blood calcium levels (hypercalcemia or hypocalcemia) brought about by a change in the secretion of PTH.

Detection of subfunctioning parathyroid glands (hypoparathyroidism) requires the use of a highly sensitive test in order to be able to measure PTH levels well below normal.^{3,4}

Hyperfunctioning of the parathyroid glands results in an increased secretion of PTH (hyperparathyroidism). Primary causes are adenomas of the parathyroid glands. In secondary hyperparathyroidism the blood calcium level is low as a result of other pathological states (e.g. vitamin D deficiency).

Today, great significance is attached to the determination of the PTH and calcium concentrations when assessing hyperparathyroidism.

The determination of PTH intraoperatively during adenoma resection in the parathyroid glands has also been reported for primary hyperparathyroidism,^{5,6,7} secondary hyperparathyroidism relating to renal failure,^{8,9} and tertiary hyperparathyroidism post renal transplant surgery.¹⁰ Because PTH has a reported half life of 3-5 minutes,¹¹ a significant drop in PTH levels after resection of the abnormal gland or glands enables the surgeon to assess the completeness of resection and whether all hyperfunctioning parathyroid tissue has been removed from the patient.¹²

The NACB guidelines recommend that baseline samples be obtained preoperation and pre-excision of the suspected hyperfunctioning gland.¹³ Specimens for PTH testing should be drawn at 5 and 10 minutes post resection and that a > 50 % reduction in PTH levels from the highest baseline be used as criteria for surgical success. Additional samples may be necessary as it has been shown that sensitivity can increase with time.¹⁴ Failure of PTH to drop below recommended levels indicates that either 1) residual hyperfunctioning tissue is still present and further exploration may be necessary, as was in the case of two patients, both with a fifth ectopic parathyroid gland requiring further surgery,⁷ or 2) a spike in PTH levels during adenoma mobilization occurred.¹⁵ Intraoperative PTH measurements

00 29/11/13

cobas

A 6229

offer fast, reliable assessment when all hyperfunctioning parathyroid tissue has been removed during the surgical process.

The Elecsys assay for determining intact PTH employs a sandwich test principle in which a biotinylated monoclonal antibody reacts with the N-terminal fragment (1-37) and a monoclonal antibody labeled with a ruthenium complex^{a)} reacts with the C-terminal fragment (38-84).

The antibodies used in this assay are reactive with epitopes in the amino acid regions 26-32 and 37-42.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: 50 µL of sample, a biotinylated monoclonal PTH-specific antibody, and monoclonal PTH-specific antibody labeled with a ruthenium complex form a sandwich complex.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode.

Reagents - working solutions

The reagent rackpack is labeled as PTH.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-PTH-Ab-biotin (gray cap), 1 bottle, 7 mL:
Biotinylated monoclonal anti-PTH antibody (mouse) 2.3 mg/L;
phosphate buffer 100 mmol/L, pH 7.0; preservative.
- R2 Anti-PTH-Ab-Ru(bpy)₃²⁺ (black cap), 1 bottle, 7 mL:
Monoclonal anti-PTH antibody (mouse) labeled with ruthenium complex 2.0 mg/L; phosphate buffer 100 mmol/L, pH 7.0;
preservative.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guideline Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

2013-07, V-23.0 English

1 / 5

PTH

Parathyroid hormone (parathormone, parathyrin) - PTH, intact

cobas[®]**Storage and stability**

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	12 weeks
on the analyzers	8 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes.

K₃-EDTA plasma.

Because of the short half-life of PTH, it is recommended that, when serum is needed, the blood should be centrifuged immediately.

Preference should be given to K₃-EDTA plasma, as it is stable longer than serum.Criterion: Method comparison serum versus plasma, slope 0.9-1.1 + intercept within $< \pm 2 \times$ analytical sensitivity (LDL) + coefficient of correlation > 0.95 .

Serum: Stable for 8 hours at 15-25 °C, 2 days at 2-8 °C, 6 months at -20 °C.

Plasma: Stable for 2 days at 15-25 °C, 3 days at 2-8 °C, 6 months at -20 °C.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 11972219122, PTH CalSet, for 4 x 1 mL
- [REF] 05618860190, PreciControl Varia, for 2 x 3 mL each of PreciControl Varia 1 and 2
- General laboratory equipment
- Elecsys 2010, MODULAR ANALYTICS E170 or **cobas e** analyzer

Accessories for Elecsys 2010 and **cobas e** 411 analyzers:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, Elecsys 2010 AssayCup, 60 x 60 reaction vessels
- [REF] 11706799001, Elecsys 2010 AssayTip, 30 x 120 pipette tips

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 03004899190, PreClean M, 5 x 600 mL detection cleaning solution
- [REF] 12102137001, AssayTip/AssayCup Combimagazine M, 48 magazines x 84 reaction vessels or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500316, Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. In exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers: PreClean M solution is necessary.

Bring the cooled reagents to approx. 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

CalibrationTraceability: This method has been standardized against a commercial PTH test (RIA). The recovery of the NIBSC 95/646 (WHO) standard was assessed by testing dilutions in human serum covering the measuring range (40-4000 pg/mL) on 16 analyzers (**cobas e** 411 and **cobas e** 601 analyzers). The mean recovery was 100 % \pm 4 %.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 7 days (when using the same reagent kit on the analyzer)
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Varia.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

MS 1137103122E260V

PTH

Parathyroid hormone (parathormone, parathyrin) - PTH, intact

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in pg/mL or pmol/L).

Conversion factors: $\text{pg/mL} \times 0.106 = \text{pmol/L}$
 $\text{pmol/L} \times 9.43 = \text{pg/mL}$

Limitations - interference

The assay is unaffected by icterus (bilirubin < 1112 $\mu\text{mol/L}$ or < 65 mg/dL), lipemia (Intralipid < 1500 mg/dL) and biotin (< 205 nmol/L or < 50 ng/mL).

The assay is affected by hemolysis ≥ 0.15 g/dL. Do not analyze samples that show visible signs of hemolysis.

Criterion: Recovery within $\pm 10\%$ of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 1500 IU/mL.

There is no high-dose hook effect at PTH concentrations up to 17000 pg/mL (1802 pmol/L).

In vitro tests were performed on 16 commonly used pharmaceuticals. No interference with the assay was found.

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

1.20-5000 pg/mL or 0.127-530 pmol/L (defined by the lower detection limit and the maximum of the master curve). Values below the lower detection limit are reported as < 1.20 pg/mL (< 0.127 pmol/L). Values above the measuring range are reported as > 5000 pg/mL (> 530 pmol/L).

Lower limits of measurement

Lower detection limit of the test

Lower detection limit: 1.20 pg/mL (0.127 pmol/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1 + 2 SD, repeatability study, n = 21).

Dilution

Not necessary due to the broad measuring range.

Expected values

15-65 pg/mL (1.6-6.9 pmol/L)^{16,17}

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents and pooled human sera in a modified protocol (EP5-A) of the CLSI (Clinical and Laboratory Standards Institute): 6 times daily for 10 days (n = 60); repeatability on MODULAR ANALYTICS E170 analyzer, n = 21. The following results were obtained:

Elecsys 2010 and cobas e 411 analyzers								
Sample	Repeatability					Intermediate precision		
	Mean		SD		CV	SD		CV
	pg/mL	pmol/L	pg/mL	pmol/L	%	pg/mL	pmol/L	%
HS ^{b)} 1	26.7	2.83	0.711	0.075	2.7	1.73	0.184	6.5
HS 2	52.5	5.56	0.853	0.091	1.6	2.07	0.220	3.9
HS 3	261	27.7	4.0	0.424	1.5	7.81	0.829	3.0

b) HS = human serum

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
Sample	Repeatability				
	Mean		SD		CV
	pg/mL	pmol/L	pg/mL	pmol/L	%
HS 1	21.9	2.32	0.44	0.05	2.0
HS 2	35.0	3.71	0.43	0.05	1.2
HS 3	123	13.04	1.31	0.14	1.1

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
Sample	Intermediate precision				
	Mean		SD		CV
	pg/mL	pmol/L	pg/mL	pmol/L	%
HS 1	23.2	2.46	0.79	0.08	3.4
HS 2	80.9	8.58	2.01	0.21	2.5
HS 3	240	25.4	6.72	0.71	2.8

Precision was determined using Elecsys reagents and controls in a separate study according to protocol EP5-A2 of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplication each for 21 days (n = 84). The following results were obtained:

Elecsys 2010 and cobas e 411 analyzers								
Sample	Repeatability					Intermediate precision		
	Mean		SD		CV	SD		CV
	pg/mL	pmol/L	pg/mL	pmol/L	%	pg/mL	pmol/L	%
PC ^{c)} Varia 1	54.6	5.79	0.657	0.070	1.2	1.11	0.118	2.0
PC Varia 2	182	19.3	2.43	0.258	1.3	3.14	0.333	1.7

c) PC = PreciControl

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers								
Sample	Repeatability					Intermediate precision		
	Mean		SD		CV	SD		CV
	pg/mL	pmol/L	pg/mL	pmol/L	%	pg/mL	pmol/L	%
PC Varia 1	54.7	5.80	0.943	0.100	1.7	0.977	0.104	1.8
PC Varia 2	184	19.5	2.65	0.281	1.4	3.11	0.330	1.7

Method comparison

A comparison of the Elecsys PTH assay (y) with a commercial PTH test (x) using clinical samples gave the following correlations (pg/mL):

Number of samples measured: 152

ms_11972103122V23.0

PTH**cobas**[®]**Parathyroid hormone (parathormone, parathyrin) - PTH, intact**Passing/Bablok¹⁸ Linear regression $y = 1.01x + 4.86$ $y = 0.83x + 19.4$ $\tau = 0.886$ $r = 0.991$

The sample concentrations were between approx. 1.4 and 1880 pg/mL (0.15 and 199 pmol/L).

Analytical specificity

No cross-reactivities were found for: Osteocalcin, PTH fragment 1-37, PTH-related protein (1-86), bone-specific alkaline phosphatase, and β -CrossLaps.

Functional sensitivity

6.0 pg/mL (0.64 pmol/L)

The functional sensitivity is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of < 20 %.

Clinical investigations in intraoperative use

In 2006, the National Academy of Clinical Biochemistry published their Laboratory Medicine Practice Guidelines for point of care testing, entitled Evidence Based Practice for Point of Care Testing.¹³ The guidelines recommend the use of intraoperative parathyroid hormone testing 1) for patients undergoing surgery for hyperparathyroidism, especially in minimally invasive or directed procedures, 2) for patients undergoing reoperation, and 3) as a replacement for traditional laboratory measurements of PTH during venous localization in order to help the angiography team guide sampling.

The guidelines further recommend for patients undergoing parathyroidectomy for hyperparathyroidism that baseline samples be obtained preoperation exploration and pre-excision of the gland, and that post-excision sampling be drawn at 5 and 10 minutes post resection with a 50 % reduction in PTH concentrations from the highest baseline level. The guidelines also caution that additional samples may be necessary.¹³

PTH testing during parathyroid surgery was conducted by several groups of investigators using the Elecsys PTH immunoassay.^{6,7,8,9,10}

The overall sensitivity and specificity of the assay to demonstrate successful surgery as defined by postoperative reduction of calcium levels was 99.6 % and 93.7 %, respectively.

References

- Silverman R, Yalow RS. Heterogeneity of parathyroid hormone: Clinical and physiologic implications. *J Clin Invest* 1973;52:1958-1971.
- Flentje D, Schmidt-Gayk H, Fischer S, et al. Intact parathyroid hormone in primary hyperparathyroidism. *Br J Surg* 1990;77:168-172.
- Nussbaum S, Potts JT. Advances in Immunoassays for Parathyroid Hormone. Clinical Applications to Skeletal Disorders of Bone and Mineral Metabolism. In Bilezikian JP, Levine MA, Marcus R (eds). *The Parathyroids: Basic and Clinical Concepts*. Raven Press, New York 1994:157-169.
- Berson SA, Yalow RS, Aurbach GD, et al. Immunoassay of bovine and human parathyroid hormone. *Proc Natl Acad Sci USA* 1963;49:613-617.
- Bergenfelz A, Nordén NE, Ahrén B. Intraoperative fall in plasma levels of intact parathyroid hormone after removal of one enlarged parathyroid gland in hyperthyroid patients. *Eur J Surg* 1991;157:109-112.
- Ohe MN, Santos RO, Kunii IS, et al. Usefulness of a rapid immunometric assay for intra-operative parathyroid hormone measurements. *Braz J Med Biol Res* 2003;36(6):715-721.
- Ohe MN, Santos RO, Kunii IS, et al. Usefulness of intra-operative PTH measurement in primary and secondary hyperparathyroidism: experience with 109 patients. *Arq Bras Endocrinol Metab* 2006;50(5):869-875.

- Seehofer D, Rayes N, Ulrich F, et al. Intra-operative measurement of intact parathyroid hormone in renal hyperparathyroidism by an inexpensive routine assay. *Langenbecks Arch Surg* 2001;386(6):440-443.
- Seehofer D, Rayes N, Klupp J, et al. Predictive value of intact parathyroid hormone measurement during surgery for renal hyperparathyroidism. *Langenbecks Arch Surg* 2005;390(3):222-229.
- Haustein SV, Mack E, Starling JR, et al. The role of intra-operative parathyroid hormone testing in patients with tertiary hyperparathyroidism after renal transplantation. *Surgery* 2005;138(6):1066-1071.
- Maier GW, Kreis ME, Renn W, et al. Parathyroid hormone after adenectomy for primary hyperparathyroidism: A study of peptide hormone elimination kinetics in humans. *Jour Clin Endocrinol Metab* 1998;83(11):3853-3856.
- Carter AB, Howanitz TJ. Intra-operative testing for parathyroid hormone: a comprehensive review of the use of the assay and the relevant literature. *Arch Pathol Lab Med* 2003;127:1424-1442.
- Nichols JH, Christenson RH, Clarke W, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines: Evidence Based Practice for Point of Care Testing. AACC Press:2006.
- Bergenfelz A, Isaksson A, Lindblom P, et al. Measurement of parathyroid hormone in patients with primary hyperparathyroidism undergoing first and reoperative surgery. *Br J Surg* 1998;85:1129-1132.
- Yang GP, Levine S, Weigel RJ. A spike in parathyroid hormone during neck exploration may cause a false-negative intra-operative assay result. *Arch Surg* 2001;136:945-949.
- Blind E. Measurement of Intact Parathyroid Hormone by an Extracting Two-Site Immunometric Assay. In: Schmidt-Gayk H, Armbruster FP, Bouillon R, (eds). *Calcium regulating hormones, vitamin D metabolites, an cyclic AMP*. Heidelberg: Springer 1990:151.
- Thomas L. Parathyroid hormone (PTH). *Clinical Laboratory Diagnosis*. TH-Books, Frankfurt. 1st english edition 1998:248-250.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT	Contents of kit
SYSTEM	Analyzers/Instruments on which reagents can be used
REAGENT	Reagent
CALIBRATOR	Calibrator
	Volume after reconstitution or mixing

COBAS, COBAS E, ELECSYS, MODULAR and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

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© 2013, Roche Diagnostics