CAT
Critically Appraised Topic

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RELEVANT EVIDENCE/REFERENCES

- 1) *Guidelines and Recommendations (most recent topics on top)*
De Pauw B et al. Revises definitions of invasive fungal disease from the European Organization for Research and Treatment of Cancer/Invasive Fungal Infections Cooperative Group and The National Institute of Allergy and Infectious Diseases Mycoses Study Group (EORTC/MSG) Consensus Group. Clin Infect Dis. 2008 June 15;46(12):1812-1821.
- 2) *Reviews*
Bourbeau PP, Ledeboer NA. Automation in clinical bacteriology. J Clin Microbiol. 2013 Jun;51(6):1658-1665.

Buchan BW, Ledeboer NA. Emerging technologies for the clinical microbiology laboratory. Clin Microbiol Rev. 2014 Oct;27(4):783-822.

Greub G, Prod'homme G. Automation in clinical bacteriology: what system to choose? Clin Microbiol Infect. 2011 May;17(5):655-660.
- 3) *Original Articles*
Baxter CG, Jones AM, Webb K, Denning DW. Homogenisation of cystic fibrosis sputum by sonication – An essential step for *Aspergillus* PCR. J Microbiol Meth. 2011 Jan 21;85:75-81.

Hammerschlag MR, Harding L, Maccone A, Smith AL, Goldmann DA. Bacteriology of sputum in cystic fibrosis: evaluation of dithiothreitol as a mucolytic agent. J Clin Microbiol. 1980 Jun;11(6):552-557.

Nebbad-Lechani B, Emirian A, Maillebuau F, Mahjoub N, Fihman V, Legrand P, Decousser JW. New procedure to reduce the time and cost of broncho-pulmonary specimen management using the Previ Isola® automated inoculation system. J Microbiol Meth. 2013 Oct 20;95:384-388.

Prattes J. Bronchoalveolar lavage fluid sample pretreatment with Sputasol® significantly reduces galactomannan levels. J Infect. 2015 May;70(5):541-552.
- 4) *Reference Works, Handbooks and Databases*
Garcia L et al. Clinical Microbiology Procedures Handbook, 3th ed. 2007. Chapter 3.2.1 Gram Stain, 3.11.2 Lower Respiratory Tract Cultures.

5) Posters, "grey literature", presentations

Bielli A, Berlingeri A, Mirone E, Landini MP. Evaluation of Copan SL-Solution for pre-treatment of lower respiratory samples for bacterial culture inoculation on the automated Walk Away Specimen Processor (WASP) [Abstract]. In: 21st European Congress of Clinical Microbiology and Infectious Diseases; 2011 May 7-10; Milan. 2011. R2467.

Boelens J, De Cooman L, Leroux-Roels I, Claeys G. Evaluation of sputum-liquefying solution for pre-treatment of respiratory samples [Abstract]. In: 23rd European Congress of Clinical Microbiology and Infectious Diseases; 2013 April 27-30; Berlin. 2013 April 27. R2752.

Boelens J, Leroux-Roels I, Claeys G. Implementation of sputum liquefying solution (SLS, Copan Italia) for use on an automated streaking robot in a routine microbiology lab [Abstract]. In: 24th European Congress of Clinical Microbiology and Infectious Diseases; 2014 May 10-13; Barcelona. 2014 May 10. R337.

Buchan B, Mackey T, Beck E, Connolly J, Ledeboer N. Evaluation of Copan SnotBuster™ for improved sputum culture results [Abstract]. In: 24th European Congress of Clinical Microbiology and Infectious Diseases; 2014 May 10-13; Barcelona. 2014 May 11. O074.

Cocquerelle V, Koebel C, Stempfer R, Graeffly C, Jehl F, Jaulhac B. Comparison of pre-treatment of sputum and endotracheal aspirate samples with SL-Solution (Copan®) or Digest-EUR (Eurobio®) [Abstract]. In: 22nd European Congress of Clinical Microbiology and Infectious Diseases; 2012 March 31-Apr 3; London. 2012 March 31. R2651.

Copan SLsolution [package insert]. Brescia: Copan Italia S.p.A. 2011.

Doyle LJ, Hall GS, Boss AC, Harrington SM. Analysis of the Copan SL Solution for liquefaction of respiratory specimens prior to Gram stain and culture [Poster]. In: 112th General Meeting of the American Society for Microbiology; 2012 June 16-19; San Francisco. 2012 June 17. C-3124.

Silbert S, Knight S, Fenton L, Widen R. Evaluation of the Copan SL Solution to improve sputum culture [Poster]. In: 112th General Meeting of the American Society for Microbiology; 2012 June 16-19; San Francisco. 2012 June 17. C-735.

Young C, Cullen K, LeBar W, Newton D. Use of COPAN SL Solution for processing sputum from patients with and without cystic fibrosis [Poster]. In: 110th General Meeting of the American Society for Microbiology; 2010 May 23-27; San Diego. 2010 May 24. C-625.

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
Copan SL solution - Product Insert & How to Use Guide

ENGLISH

INTENDED USE

Copan SLsolution (Sputum Liquefying solution) is a ready-to-use mucolytic agent for liquefying sputum specimens prior to planting and streaking for the isolation of bacteria and fungi that cause respiratory tract infections..

SUMMARY

Copan SLsolution consists of phosphate buffered solution of dithiothreitol (DTT). The DTT solution is provided in a ready-to-use liquid format and the intended use is for digesting sputum samples prior to culturing. Copan SLsolution has been validated for use with Copan Walk-Away Specimen Processor (WASP) an automated system for planting and streaking of specimens onto culture media plates..

Copan SLsolution is available in different configurations: in bulk packs comprising multiple screw cap tubes filled with 1 ml of DTT solution or as kits comprising of one screw cap tube containing 1 ml of DTT solution plus a sterile disposable Pasteur pipet or Sputum Dipper transfer device. All SLsolution configurations are packaged within barrier plastic film pouches that have a modified atmosphere to maintain the stability and activity of the DTT reagent during the products shelf life.

REAGENTS

SLsolution components:

DTT (dithiothreitol) 1mg
PBS (phosphate buffer solution) 1ml

PRECAUTIONS

- For *in vitro* diagnostic use.
- Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
- All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimen, containers and media after their use.
- Directions should be read and followed carefully.

STORAGE

This product is ready to use and no further preparation is necessary. The unopened bulk pack or individual kit pack can be stored at 5 - 25°C until used or until the expiration date. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box. The packaging for SLsolution comprises of barrier plastic film together with a modified atmosphere inside package maintains the stability of the product. After opening the package the tube(s) of SLsolution should be used immediately or must be stored at 2-8°C and used within 24 hours. Slight yellow coloration of the DTT reagent does not indicate any product deterioration.

PRODUCT DETERIORATION

Do not use **Copan SLsolution** if: (1) the product shows visible marks of damage or contamination; (2) there is evidence of leakage; (3) the expiration date has passed; (4) the package is open or more than 24hr and has not been stored in a refrigerator (5) there are other signs of deterioration.

MATERIALS SUPPLIED

Catalog No.	Product Descriptions	Pack Size	Suitable for Automation
097CE 097CE.A	Copan SLsolution in bulk: - 12X80 mm polypropylene screw-cap tubes with internal conical shape filled with 1ml of DTT solution in liquid phase	30 units per shelf pack 10x30 units per box	YES
095CE 095CE.A	Copan SLsolution Kit: - One 12X80 mm polypropylene screw-cap tube with internal conical shape filled with 1ml of DTT solution in liquid phase - One sterile pasteur pipet for sputum sample transfer	20 units per shelf pack 10x20 units per box	YES
099CE.A	Copan SLsolution in bulk: - 12X80 mm polypropylene screw-cap tubes with internal conical shape filled with 0.5 ml of DTT solution in liquid phase	30 units per shelf pack 10x30 units per box	YES

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for the cultivation and isolation of bacteria and fungi. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques

INSTRUCTIONS FOR USE

SLsolution in bulk:

- Open the plastic pack containing the tubes of SLsolution remove one tube and unscrew the cap
- Transfer sputum sample into the opened tube. NOTE: best results can be achieved transferring the sample into the SL solution respecting volumetric ratio about 1:1.
- Re-cap the tube of SLsolution
- Leave the test tube at room temperature for at least 15 minutes. A prolonged contact time, up to 6 hours does not impair the survival of bacteria or fungi contained in the sample. Vortex the test tube for 30 seconds at 2000/2500 rpm or until the mixture has liquefied.
- NOTE: for especially viscous samples, a minimum contact time above 15 minutes may be necessary. The operator is to ensure that the sample has actually liquefied
- Aseptically remove aliquots of the liquefied sputum and inoculate onto appropriate bacteriology culture medium.

SLsolution Kit with sterile Pasteur pipet:

- Peel open one sterile graduated pipet
- Suck up about 1ml of sputum specimen into the pipet. Open the plastic pack containing the tube of SLsolution remove one tube and unscrew the cap
- Transfer sputum sample into the opened tube.
- Re-cap the tube of SLsolution
- Leave the test tube at room temperature for at least 15 minutes. A prolonged contact time, up to 6 hours does not impair the survival of bacteria or fungi contained in the sample. Vortex the test tube for 30 seconds at 2000/2500 rpm or until the mixture has liquefied.
- NOTE: for especially viscous samples, a minimum contact time above 15 minutes may be necessary. The operator is to ensure that the sample has actually liquefied..
- Aseptically remove aliquots of the liquefied sputum and inoculate onto appropriate bacteriology culture medium.

LIMITATIONS

- In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens.
- Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.
- Copan SLsolution is not tested for acid resistant bacilli culture.
- Performance testing with Copan SL solution was conducted using laboratory strains spiked into the SL Solution tube and not using human specimens.
- Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals. Specimens should be collected as soon as possible after the clinical onset of disease. Highest bacterial titers are present during the acute illness.

WARNINGS

- Do not re-pack.
- Not suitable for any other application than intended use.
- The use of this product in association with any diagnostic assay or with any diagnostic instrumentation should be validated by the user before using.
- Do not use if the product is visibly damaged
- Do not ingest the medium.
- Directions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
- Product to be handled by trained personnel only.
- It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes must be disposed of according to laboratory regulations for infectious waste.
-  Copan SL solution is for single use only; reuse may cause a risk of infection and/or inaccurate results.
- Copan SL solution is for in-vitro diagnostics use only and is in no way intended for a curative or prophylactic purposes.

WASTE DISPOSAL

Unused reagents may be considered as non hazardous waste and disposed of accordingly.

Dispose of used reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

PERFORMANCE CHARACTERISTICS:

An investigation into the survival of respiratory pathogens listed below in contrived specimens that have been stored for 6 hours at room temperature following homogenisation with SL solution showed that organism remained viable and reculturable:

Pseudomonas aeruginosa ATCC 27853
Haemophilus influenzae ATCC 10211
Stafilococcus aureus ATCC 6538

Streptococcus pneumoniae ATCC 6305
Streptococcus pyogenes ATCC 19615
Moraxella catarrhalis ATCC 25238
Candida albicans ATCC 10231

PERFORMANCE TEST RESULTS

STRAIN	CFU ZERO TIME COUNT	CFU COUNT AFTER 6 HOURS AT ROOM TEMPERATURE
Pseudomonas aeruginosa ATCC 27853	396	167
Haemophilus influenzae ATCC 10211	320	88
Staphylococcus aureus ATCC 6538	285	270
Streptococcus pneumoniae ATCC 6305	350	184
Streptococcus pyogenes ATCC 19615	466	178
Moraxella catarrhalis ATCC 25238	501	37
Candida albicans ATCC 10231	275	260

Procedure for testing SL solution:

- Starting from a fresh culture plate of tester organism prepare 0.5 Mc Farland inoculum of the tester organism in PBS
- prepare appropriate dilutions in PBS of original 0.5 McF inoculum in order to obtain a diluted inoculum containing from 2000 to 20000 CFU/100 ul
- with a micropipette spike 100 ul of chosen dilution into SLsolution tube
- vortex inoculated tube for 30 seconds
- spread 100 ul of inoculated SL solution onto appropriate culture plate
- incubate the plate at 35°C for 48 hours in order to obtain zero time count plate
- hold the inoculated SL tubes for 6 hours at room temperature
- vortex inoculated tube for 30 seconds
- spread 100 ul of inoculated SL solution onto appropriate culture plate
- incubate the plate at 35°C for 48 hours in order to obtain six hours time count plate
- vortex for 30 seconds and spread 100 ul to obtain six hours plate and incubate it at 35°C for 48 hours

ACCEPTABILITY LIMITS : >10 CFU/plate after six hours at room temperature.

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BIBLIOGRAPHY

- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger and W.C. Winn, Jr. 1992. Color Atlas and Textbook of Diagnostic Microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality Control of Microbiological Transport Systems. Approved Standard M40-A.
- 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Isenberg, H.D., 1998. Essential Procedures for Clinical Microbiology. Chapter 14.12, Page 787. Packaging and Shipping Infectious Substances.
- National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents. Approved Standard H5-A3.

Tabella dei Simboli/Index of Symbols/Table des Symboles/Symbole/Tabla de Símbolos/ Quadro de Símbolos

Symbol/Symbole/Simbolo	Significato/Meaning Signification/Bedeutung/Significado/ Significado
	Fabbricante/ Manufacturer/Fabricant/Hersteller/Fabricante/ Fabricante
	In Vitro Diagnostic
	Non riutilizzare/Do not reuse/Ne pas réutiliser/Nicht zur Wiederverwendung/No reutilizar/ Não voltar a usar
	Numero di catalogo/Catalogue number/Référence du catalogue/Bestellnummer/Número de catálogo/ Referência do catálogo
	Limiti di temperatura/Temperature limitation/Limites de temperature/Temperaturbegrenzung/Limites de temperatura/ Limites de temperatura
	Utilizzare entro/Use by/Utiliser jusque/Verwendbar bis/Fecha de caducidad/ Prazo de validade
	Consultare le istruzioni per l'uso/Consult Instructions for Use/Consulter les instructions d'utilisation/Gebrauchsanweisung beachten/Consulte las instrucciones de uso/ Consultar as instruções de utilização
	Strappare per aprire/Peel/Décoller/Abziehen/Desprender/ Destacável
	Codice del lotto (partita)/ Batch code (Lot)/Code de lot (Lot)/Chargencode (Chagenbezeichnung)/Código de lote (Lote)/ Código do lote (Lote)
	Contenuto sufficiente per <n> test/Contains sufficient for <n> tests/Contenu suffisant pour <n> tests/Ausreichend für <n> Tests/Contenido suficiente para <n> pruebas/ Contém suficiente para <n> testes
 Copen Italia Spa Via Perotti, 10 Brescia, Italy	 Copen Italia Spa Via Perotti 10 25125 Brescia Italy Tel: +39 030 2687211 Fax: +39 030 2687250 E-mail: info@copenitalia.com Website: www.copeninnovation.com